

ETO - 4.5 (124 LITERS) FULLY AUTOMATIC OPERATED



NEMR Industries Pvt. Ltd.
Empowering Nation with Quality Products



TECHNICAL SPECIFICATION

We are an eminent Manufacturer and Supplier of ETO Sterilizers from Himachal Pradesh, India. The ETO Sterilizers are user-friendly, compact in design low priced, and easy maintenance. We have become the primary choice of the buyers as we deliver even the bulk orders on time and all the credit goes to our logistics.

Particular uses are sectors such as pharmaceutical, food plant, healthcare, Hospital & Medical College and other industries.

MATERIAL OF CONSTRUCTION:

Inner Chamber S.S 316. Water

Jacket S.S 304 Vacuum Pump:

Vane Type Oil Type

Control: Individual push-button for vacuum pump, gas, Puncturing, vacuum breaker, and heating with digital PID controller

Operating voltage: 230 voltage /50 Hz

Non-Hazardous: Total cycle under negative Pressure

TECHNICAL SPECIFICATION:

Rectangular chamber, water, pipeline made of SS 304 quality

Working Temperature: 37°C to 55°C Mode of

Operation: Manual/Automatic Gas Injection: Cartridge

Insulation: Glass Wool 25 - 50mm Full

Inner chamber SS 316/304

Enclosure: Cabinet Type with Powder Coated

Temperature Controller: Digital types LED Display. We are an eminent Manufacturer and Supplier of ETO Sterilizers from Himachal Pradesh, India. The ETO Sterilizers are user-friendly, compact in design and low priced. Temperature sensor RTD-100 Type and easy maintenance. We have become the primary choice of the buyers as we deliver even the bulk orders on time and all the credit goes to our logistics. Particular Door: Single Point Locking System/Radian arrangement locking system. uses are sectors such as pharmaceutical, food plant, healthcare, Hospital & Medical College and other industries.

OPTIONAL SPARE ACCESSORIES

- Vacuum gauge
- Temperature Sensor
- ETO Cartage 100ml & as per size
- Joint less Gaskets: - silicon food grade.
- biological indicator, packing roll, sealing machine
- Sealing /packing machine
- Packing Roll all sizes.

AVAILABLE SIZE

Model	Inner Size mm	Volume ltr	Power Kw	Voltage
NEP-ETO-509	305X305X750	70	1.0	230V
NEP-ETO-509	305X450X610	96	1.5	230V
NEP-ETO-509	450X450X610	124	2.0	230V
NEP-ETO-509	450X450X915	185	2.5	230V
NEP-ETO-509	450X450X1220	247	3.0	230V
NEP-ETO-509	610X610X915	340	3.0	230V
NEP-ETO-509	610X610X1220	454	4.0	230V