

Fume Hood has heavy Stainless Steel Construction 304 Grade Size: 4'x2'x2', 6X2X2, Work table is made of SS-316 or ceramic glazed tiles, a small S.S. Sink, Water Tap, Water outlet provided. A sliding glass door Moves vertically with counter weight. Complete with fluorescent light. Centrifugal motor fan with S.S. Fitted with the exhaust fumes will vented with PVC pipe outlet. (160 diameter).

National Fume Hood has heavy Stainless Steel Construction 304 Grade Size: 4'x2'x2', 6X2X2, Work table is made of SS-316 or ceramic glazed tiles, a small S.S. Sink, Water Tap, Water outlet provided. A sliding glass door Moves vertically with counter weight. Complete with fluorescent light. Centrifugal motor fan with S.S. Fitted with the exhaust fumes will vented with PVC pipe outlet. (160 diameter)

Specifications:

Stainless Steel Construction of 304 grade Available in size of 2' x 2' Work table comprising ceramic glazed tiles Provision of stainless steel sink, water tap and outlet Sliding glass door that provides vertical movements with counter weight.

Fluorescent light and centrifugal motor fan support, SS Sink, SS tab, etc. Unit is widely demanded by the pharmaceutical industry as well as in huge chemical laboratories.

- ❖ Compact design
- ❖ Low noise Level
- ❖ Excellent air Exit flow through duct



Overall Dimensions / Working Size

- ❖ The system should have following Overall Dimensions:
- ❖ Overall Length of Fume Hood: 1 500 - 1 550 mm
- ❖ Overall Width of Fume Hood: 750 - 1500 mm
- ❖ Overall Height of Fume Hood: 1500 - 2500 mm
- ❖ Length of Base Cabinet: 1000 - 1500 mm
- ❖ Height of Base Cabinet: 700 - 800 mm.

Body Features.

Double Wall Construction o Body thickness: 10 mm (Min.) o Completely made from GI sheet with Highly corrosion resistant epoxy powder coating r Inner Chamber - Chemical & Heat Resistance, Fire retardant, smooth finish, easily cleanable, made out of durable PRL sheets of thickness 5 mm (Min.)/ SS 304 of 18- 20 gauge thickness. r Should be provided with Fume Hood installation Kit and Accessories we provided with Safety Device Trip. Working Table Top Granite / M.S Powder Coated Sheet Covered with P.P Sheet/ ss 304 Thickness of granite 18 mm (Min.).

IQ Q/OQ On site IQ, OQ of instrument along with document to be provided & supplier to assist till satisfactory PQ of instrument.

Door / Sash/ Shutter Thickness -4 mm (min.) Material - Toughened Glass Door vertical Folding Type with adjustable height